

ADVANTAGES OF THE MOOSE

- ▶ After familiarisation patient can confidently mobilise solo, and more often, with all necessary equipment
- ▶ Physical therapist can concentrate on teaching exercises and assessing patients
- ▶ Nurses can concentrate on nursing and caring for patients holistically
- ▶ Frees up hospital beds earlier, shortening waiting lists
- ▶ Reduces preparation time for mobilising patients from 7 min to 1 min
- ▶ Saves room space. All equipment is reduced to one unit
- ▶ Stackable. Saves storeroom space when not in use
- ▶ Height adjustment ranges from small child to tall adult. Suits all ages.

INDICATIONS FOR USE

- ▶ All major surgical patients
- ▶ Patients requiring 2 conventional IV Poles
- ▶ Patients requiring walking frame plus IV or drains
- ▶ Patient requiring oxygen and other equipment while mobilising
- ▶ Any patient having stability problems with present outdated equipment

HOME USE

Patients who previously could not be discharged because they required too much separate equipment e.g. oxygen and walking frame.

COMING SOON

Dedicated home-use model will be available soon.

AboutTime
TECHNOLOGIES

les@abouttimetechnologies.com
0400 273 789
russell@abouttimetechnologies.com
0438 878 424

71 Goongarrie Drive
COOLOONGUP 6168
T: +61 8 9592 4673
F: +61 8 9529 3004

www.abouttimetechnologies.com

Its About Time.



The Patient Mobilizer has arrived.

moose
mobilizer 



Invented by a registered nurse with a surgical background.

This patient mobiliser is used to mobilise major surgical patients sooner and more frequently, giving them a shorter recovery period.

Four pieces of hospital equipment:

- ▶ 2 IV poles
- ▶ 1 walking frame
- ▶ 1 oxygen holder
- ▶ plus 2 baskets for drains and accessories


All superior to the original equipment.



INTEGRATED ELECTRICS
5 outlets for 5 pumps.
One cord required to wall outlet.

STACKABLE
When not in use frees up storage space.

VERY STABLE BASE
Manoeuvrable independent wheels free walking area.
No castors in the way.



DETACHABLE EXTENSIONS
Remove when not required for gravity feed infusions.



FIXED HEIGHT INFUSION HOOKS
Most infusions should be on these hooks. Some pumps may require gravity assist. Use extension.



ADJUSTABLE INTERCHANGABLE HANDLES
Children to tall adults in one easy adjustment.

ANGLED MAIN PUMP AREA
No bending to read or operate.



O2 CYLINDER ON BOARD
Also helps lower centre of gravity when in situ.

VERY LOW CENTRE OF GRAVITY
all pumps, equipment and fluids lowered.